

# Desktop Fast Autoclave Sterilizer 35L / 50L / 90L

Item Number: KT-B04



## Introduction

The desktop fast steam sterilizer is a compact and reliable device used for rapid sterilization of medical, pharmaceutical, and research items. It efficiently sterilizes surgical instruments, glassware, medicines, and resistant materials, making it suitable for various applications.

[Learn More](#)

Model	KT-XD35DV	KT-XD50DV	KT-XD90DV
Sterilization chamber volume	35L $\phi$ 300×500mm	50L $\phi$ 340×550mm	90L $\phi$ 400×720mm
Fixed work stress	0.22Mpa		
Rated working temperature	134°C		
Sterilization temperature selection	105-134°C		
Timing time selection range	0-99min		
Heat uniformity of sterilization room	$\leq \pm 1^\circ\text{C}$		
Power/supply voltage	AC220V 50Hz/2.7KW		AC220V 50Hz/4KW
Overall dimensions (mm)	730×550×510	770×600×540	1000×690×610
Packing size (mm)	810×680×650	880×730×700	1070×740×700
Gross weight/net weight	115/93Kg	130/106Kg	200/180Kg