

Four-Body Horizontal Jar Mill

Item Number: KT-HJM



Introduction

The four-body horizontal tank mill ball mill can be used with four horizontal ball mill tanks with a volume of 3000ml. It is mostly used for mixing and grinding laboratory samples.

[Learn More](#)

Processing principle	Friction
Application sample features	Fine Medium to low hardness Brittle Dry or low viscous
Processing type	grinding mixing
Maximum injection size	< 5 mm
Sample particle size range	1-20um
Number of ball mill jars	4
The largest ball milling tank volume available	3000ml
Available minimum milling jar volume	250ml
Ball mill tank optional material	Zirconia Tungsten Carbide Agate Nylon Stainless Steel Polyurethane etc
Grinding media optional type	Spherical Cylindrical
Grinding ball optional material	Zirconia Tungsten Carbide Agate Stainless Steel Corundum
Grinding method	Dry grinding wet grinding
Number of rubber rollers:	3 sticks
Rubber roller speed	900 r/min
Rubber roller cover material	Polyurethane
Rubber roller movement method	electric control
net weight	88kg
Dimensions	670*660*310mm
Protection class	IP65
standard	CE