

## Horizontal Autoclave Steam Sterilizer (Microcomputer)

Item Number: KT-C06



### Introduction

The horizontal autoclave steam sterilizer adopts the method of gravity displacement to remove the cold air in the inner chamber, so that the content of steam cold air in the inner chamber is less, and the sterilization is more reliable.

[Learn More](#)

Technical data	KTM -150YDB	KTM -200YDB	KTM -280YDB
Sterilizing chamber volume	150L φ440×1000mm	200L φ515×1000mm	280L φ600×1000mm
working pressure	0.22 MPa		
working temperature	134 degree		
Adjustment of temperature	40 - 134 degree		
Time of sterilization	0-99 min		
Time of drying	0-99 min		
Heat average	≤± 2 degree		
Power	9KW/380V 50Hz	9KW/380V 50Hz	12KW/380V 50Hz
Dimension	1400×650×1600	1400×750×1700	1520×910×1900
Transport dimension	1560×820×1820	1560×910×1880	1680×1080×2100
G.W/N.W	430/340 KG	436/350 KG	570/462 KG